



## Diode Laser Permanent Hair Removal Machine 755nm 808nm 1064nm Medical Beauty Equipment

### Our Product Introduction

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: TM-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



#### Product Specification

- Name: Diode Laser Permanent Hair Removal Machine
- Machine App: Permanent Laser Hair Removal Machines
- Machine Power: 1200W
- Coolingtype: TEC+Sapphire Contact Cooling
- Laser Type: Other
- Style: Stationary, Vertical
- Energy Density: Up To 120 J/cm2
- Display: 15-inch Touch Screen
- Pulse Duration: 10-400 Ms
- Highlight: Diode Laser Permanent Hair Removal Machine 755nm  
,  
Diode Laser Permanent Hair Removal Machine

#### Instrument Product Component

The internal hardware assembly and design are the first to join the aluminum alloy material, which can greatly extend the life of the instrument and play the function of the hardware itself, and can better protect the core components from accidental damage.



## Product Description

### Permanent Diode Laser Hair Removal Machine 755 808 1064Nm Medical Beauty Equipment

#### Product Description:

Designed to cater to all skin types ranging from I-VI, the Hair Removal Laser Machine is guaranteed to give you the results you desire without causing any harm to your skin. The laser treatment is gentle yet effective, making it a popular choice among individuals who want to get rid of unwanted hair without experiencing any pain or discomfort.

The Hair Removal Laser Machine is equipped with several safety features such as Skin Contact Cooling and Automatic Power Cut-off. These features ensure that your skin remains safe during the entire treatment process. The Skin Contact Cooling feature helps to reduce the discomfort that may be caused by the laser treatment, while the Automatic Power Cut-off ensures that the machine shuts off automatically in case of any malfunction.

Using the Hair Removal Laser Machine is easy and convenient. The machine is designed to provide fast and efficient results, making it perfect for individuals who are always on-the-go. With regular use, you can expect to see a significant reduction in hair growth, leading to smoother and more beautiful skin.

Whether you are tired of shaving or waxing, the Hair Removal Laser Machine is the perfect solution for you. It is an affordable and effective alternative to traditional hair removal methods, making it a popular choice among individuals who want to save time and money.

In conclusion, if you want to achieve smooth and hair-free skin, the Hair Removal Laser Machine is the perfect choice for you. It is a Smooth Skin Laser Machine that utilizes a powerful 600W Laser Type to provide fast and efficient results. With its safety features and easy-to-use design, you can rest assured that your skin will remain safe and beautiful throughout the entire treatment process.

#### Features:

Product Name: Hair Removal Laser Machine

Style: Stationary, Vertical

Machine App: Permanent Laser Hair Removal Machines

Coolingtype: TEC+Sapphire Contact Cooling

Energy Density: Up To 120 J/cm<sup>2</sup>


Safety Features:

Skin Contact Cooling

Automatic Power Cut-off

#### Technical Parameters:

Hair Removal Laser Machine	
Hair Depilator Machine	
Painless Hair Removal System	
Cooling Fan	2pcs Fans
Machine Power	1200W
Style	Stationary, Vertical
Power	600W
Pulse Duration	10-400 Ms
Machine App	Permanent Laser Hair Removal Machines
Skin Types	I-VI
Laser Type	Other
Safety Features	Skin Contact Cooling, Automatic Power Cut-off
Display	15-inch Touch Screen



## More than just display

2.4-inch high-definition LED touch screen,  
treatment parameters can be adjusted directly  
on the handle screen.



- ✓ Energy adjustment
- ✓ Frequency adjustment
- ✓ Pulse width adjustment
- ✓ Real-time temperature monitoring



Ice Hair Removal



Selectable  
Spot Size



Diode Laser



Fast Light



Adjustable  
Frequency



Adjustable Energy

01



### LED touch screen

LED display, large screen  
operation and sensitive touch

02

### Vent design

Back air cooling + full radiator  
with multiple vents for more  
peace of mind



## Instrument Product Component

The internal hardware assembly and design are the first to join the aluminum alloy material, which can greatly extend the life of the instrument and play the function of the hardware itself, and can better protect the core components from accidental damage.





# PRODUCT DETAILS

Laser type	Diode laser
Laser wavelength	808nm/755nm/1064nm /940nm
Pulse frequency	1~10Hz
Spot size	12*12mm/12*35mm
Pulse Width	1-400 milliseconds (adjustable)
skin type	I-VI skin types
Laser module power	600W/800W/1000W11200W/1600W
Voltage 电压	110V/220V+20V 50/60HZ
Cooling system	Air + water + compressor cooling
Energy	1-100j

## Applications:

The LZR TM-1 Hair Removal Laser Machine is designed for use in various settings and occasions. It can be used in salons and spas, as well as at home. This hair removal laser machine is perfect for people who want to save time and money on frequent salon visits. With its powerful 600W motor and energy density of up to 120 J/cm<sup>2</sup>, it is capable of removing hair from all skin types, from I to VI.

The LZR TM-1 Hair Removal Laser Machine is also equipped with a TEC+Sapphire contact cooling system that ensures comfortable and pain-free hair removal. Whether you are removing hair from your legs, arms, underarms, or bikini area, you will experience a painless hair removal system that keeps your skin cool and comfortable.

This hair removal laser machine is designed for stationary and vertical use. Its stylish and modern design will complement any salon or spa décor. It is also lightweight and easy to move around, making it perfect for use at home or on the go.

The LZR TM-1 Hair Removal Laser Machine is perfect for use on all occasions and scenarios. Whether you are getting ready for a special event, or just want to feel confident in your own skin, this hair depilator machine has got you covered. With its advanced features, painless hair removal system, and versatile design, you can achieve silky smooth skin anytime, anywhere.

## Customization:

LZR TM-1 Hair Removal Laser Machine is a powerful and effective permanent laser hair removal machine that can remove unwanted hair from all skin types (I-VI). Its energy density of up to 120 J/cm<sup>2</sup> and power of 600W make it an ideal hair removal solution. The laser type used in this machine is other and it is made in China Guangzhou.

Our product customization services allow you to tailor the LZR TM-1 Hair Removal Laser Machine to your specific needs. You can choose from a range of options to create your own painless hair removal system. We can help you add features such as adjustable energy density, customizable laser types, and more.

Our unwanted hair remover device is designed to provide you with a comfortable hair removal experience. With our painless hair removal system, you can say goodbye to the pain and discomfort associated with traditional hair removal methods.

Trust us to provide you with the best product customization services for your LZR TM-1 Hair Removal Laser Machine. Contact us today to learn more about our customization options.

## Support and Services:

The Hair Removal Laser Machine is a state-of-the-art device designed to provide effective hair removal treatment at home. Our product technical support and services include:

Online user manual and instructional videos to help you set up and use the machine

24/7 customer support for any technical or operational issues you may encounter

A one-year warranty covering any defects or malfunctions in the machine

Regular software updates to ensure the machine is always up-to-date and functioning optimally

Access to our online community where you can connect with other users and share tips and advice

We are committed to providing our customers with the best possible experience using our Hair Removal Laser Machine. If you have any questions or concerns, please do not hesitate to contact us.

## FAQ:

Q: What is the brand name of this hair removal laser machine?

A: The brand name is LZR.

Q: What is the model number of this hair removal laser machine?

A: The model number is TM-1.

Q: Where is this hair removal laser machine manufactured?

A: It is manufactured in Guangzhou, China.

Q: What type of hair can this laser machine remove?

A: This laser machine can remove hair of all colors and types, including dark and coarse hair.

Q: Is this laser machine safe to use on all skin types?

A: Yes, this laser machine is safe to use on all skin types, but we recommend consulting with a dermatologist before using it on sensitive skin.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District,  
Guangzhou