



## Customized 800W Diode Hair Removal Laser Machine 60 Million Light Outputs

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: TM-1
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



#### Product Specification

- Power: 600W
- Name: Diode Hair Removal Laser Machine
- Skin Types: I-VI
- Energy Density: Up To 120 J/cm<sup>2</sup>
- Coolingtype: TEC+Sapphire Contact Cooling
- Style: Stationary, Vertical
- Safety Features: Skin Contact Cooling, Automatic Power Cut-off
- Cooling Fan: 2pcs Fans
- Laser Type: Other
- Highlight: Customized Diode Hair Removal Laser Machine, 800W Diode Hair Removal Laser Machine



## Product Description

### Customizable 800W Diode Laser Hair Removal Machine 60 Million Light Outputs

#### Product Description:

The Hair Depilator Machine is designed to cater to a wide range of skin types. It is suitable for skin types I-VI, making it a versatile option for people with different skin tones. It offers a pulse duration of 10-400 Ms, which allows it to target hair follicles with precision and effectiveness. This laser hair removal device is a highly effective solution for people who are looking for a permanent hair removal device. Whether you're looking to remove hair from your legs, arms, underarms, or bikini line, the Hair Depilator Machine is the perfect device for you. It is easy to use and provides long-lasting results. The machine uses laser technology to target and destroy hair follicles, preventing hair growth in the treated areas. With consistent use, this Permanent Hair Removal Device can help you achieve smooth and hairless skin.

#### Features:

Product Name: Hair Removal Laser Machine

Style: Stationary, Vertical

Safety Features: Skin Contact Cooling, Automatic Power Cut-off

Display: 15-inch Touch Screen

Cooling Fan: 2pcs Fans

Power: 600W

This Laser Hair Treatment Device is a professional hair removal device that comes with safety features like skin contact cooling and automatic power cut-off. With a 15-inch touch screen display, it is easy to use and control. The 2pcs fans provide cooling for the device while it is in use. The device has a power of 600W and is designed to be used in a stationary, vertical style.



#### Technical Parameters:

Style:	Stationary, Vertical
Skin Types:	I-VI
Machine App:	Permanent Laser Hair Removal Machines
Display:	15-inch Touch Screen
Power:	600W

Pulse Duration:	10-400 Ms
Cooling Type:	TEC+Sapphire Contact Cooling
Safety Features:	Skin Contact Cooling, Automatic Power Cut-off
Cooling Fan:	2pcs Fans
Machine Power:	1200W



## COMPARISON BEFORE AND AFTER USE





### no-invasive 1064

(Washing red eyebrows/colorful tattoos, etc.)

Mainly targeting warm color pigments



### no-invasive 532

(Wash blue eyebrows/tattoos/eyeliner)

Mainly targeting cold color schemes and black pigments



### laser 532

(Deep eyebrow washing/  
tattoo washing/mole spotting)

Mainly targeting cold color schemes and black pigments



### laser 1064

(Deep brow wash spots/birthmarks/tattoos)

Mainly targeting warm color pigments



### laser 1320/Black faced doll

(Oil control/Pore cleaning/Blackhead removal/  
Carbon powder application required)

Mainly (cleaning superficial pigments  
in the epidermis)

## Applications:

The LZR Hair Removal Laser Machine is suitable for both home and salon use. It can be used to remove hair from various body parts, including the face, arms, legs, underarms, and bikini area. The device is ideal for individuals with busy schedules who want to remove hair quickly and efficiently. The machine is also perfect for salon owners who want to offer laser hair treatment services to their clients.

The LZR Hair Removal Laser Machine features TEC+Sapphire Contact Cooling, which ensures that the skin remains cool during the hair removal process. The cooling system helps to prevent skin damage, making the device suitable for all skin types. The machine has a power output of 1200W, making it efficient and fast.

The LZR Hair Removal Laser Machine is ideal for individuals who want to achieve permanent hair reduction without experiencing any pain or discomfort. The machine is also perfect for individuals who are tired of using traditional hair removal methods, such as waxing and shaving. With the LZR Hair Removal Laser Machine, you can achieve smooth, hairless skin in a matter of minutes.

In summary, the LZR Hair Removal Laser Machine is a reliable and effective hair removal device that is suitable for both home and salon use. It is a painless hair removal system that utilizes advanced laser technology to offer permanent hair reduction. The machine features TEC+Sapphire Contact Cooling, a 15-inch touch screen display, and a power output of 1200W, making it efficient and fast.

## Customization:









### Support and Services:

The Hair Removal Laser Machine is a high-tech device that is designed to provide a painless and effective way to remove unwanted hair from various parts of the body. Our product technical support and services include:

Online user manual and instructional videos to help you get started with your device

Expert technical support to answer any questions or concerns you may have about your device

Regular software updates to ensure optimal performance and safety

Repair services in case of any issues with the device

We are committed to providing our customers with the highest quality of technical support and services to ensure that you get the most out of your Hair Removal Laser Machine.

### FAQ:

Q: What is the Brand Name of this Hair Removal Laser Machine?

A: The Brand Name of this Hair Removal Laser Machine is LZR.

Q: What is the Model Number of this Hair Removal Laser Machine?

A: The Model Number of this Hair Removal Laser Machine is TM-1.

Q: Where is this Hair Removal Laser Machine manufactured?

A: This Hair Removal Laser Machine is manufactured in China Guangzhou.

Q: What is the wavelength of the laser used in this Hair Removal Laser Machine?

A: The wavelength of the laser used in this Hair Removal Laser Machine is 808nm.

Q: What is the expected lifespan of the laser in this Hair Removal Laser Machine?

A: The expected lifespan of the laser in this Hair Removal Laser Machine is 10 million shots.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou



